

HUGO LP-HLED-C:

'HUGO' Ceiling Mounted LED Examination Lamp

www.lipmedical.com



The future looks good.


LIP MEDICAL
HEALTHCARE DESIGN

Specifications:

PR-HLED-C

Lamp Type	5 xLED module, 12W
Light Intensity	50,000 Lux@0.5m
Light Field Diameter	Ø150mm@0.5m
Colour temperature	4,500 Kelvin
LED Life Expectancy	>45,000 hours
Control	Illuminated On/Off Switch mounted on the lamp head
Safety Standards	EN60601-1, EN60601-2-41
Arm Length	1540mm
Supply Voltage	230V, 50-60Hz
Total Weight	Approx. 7.5kg
Class of Protection	Class II
Floor to Ceiling Height	
Range:	2500mm - 2700mm
	2700mm - 2900mm
	2900mm - 3100mm



T: 0113 251 5119
F: 0113 251 5100
enquiries@lipmedical.com

Offices

1200 Century Way, Thorpe Park Business Park, Colton, Leeds, LS15 8ZA
Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA

