

**HUGO LP-HLED:**

# 'HUGO' Wall Mounted LED Examination Lamp

www.lipmedical.com



The future looks good.

**LIP MEDICAL**  
HEALTHCARE DESIGN

## Specifications:

# PR-HLEDWR – Wall Mounted LED Examination Light

<b>Lamp Type</b>	5 xLED module, 12W
<b>Light Intensity</b>	50,000 Lux@0.5m
<b>Light Field Diameter</b>	Ø150mm@0.5m
<b>Colour temperature</b>	4,500 Kelvin
<b>LED Life Expectancy</b>	>45,000 hours
<b>Control</b>	Illuminated On/Off Switch mounted on the lamp head
<b>Safety Standards</b>	EN60601-1, EN60601-2-41
<b>Arm Length</b>	1200mm
<b>Supply Voltage</b>	230V, 50-60Hz
<b>Total Weight</b>	Approx. 2kg
<b>Class of Protection</b>	Class II



**T:** 0113 251 5119  
**F:** 0113 251 5100  
[enquiries@lipmedical.com](mailto:enquiries@lipmedical.com)

### Offices

1200 Century Way, Thorpe Park Business Park, Colton, Leeds, LS15 8ZA  
Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA

