

Mediline, Horizontal Bedhead Units



MANUFACTURED IN THE UK


LIP MEDICAL
HEALTHCARE DESIGN



Professional Healthcare Design.

www.lipmedical.com



T: 01244 288 299
enquiries@lipmedical.com

Our primary focus is the understanding of and commitment to enhance the patient environment.

With a broad range of knowledge and experience, we feel passionately about both design and functionality of the patient environment.

**We are
LiP Medical**

Healthcare Design.

Mediline



Supply Units

Horizontal Bedhead Supply Unit

The Mediline horizontal bedhead unit offers an ergonomic, compact stylish solution to deliver essential services to the patient. The units are configured to suit each individual installation or project requirement with modules supplied fully assembled as standard. Electrical sockets, nurse call, data and medical gases are fitted to meet each individual client specification.

Electrical sockets are fitted to individual accessory plates and can be fully wired and tested with essential/non-essential sockets available for standard ward areas and blue isolated earth sockets for A&E, critical care or operating departments.

Medical gas terminal units from all manufacturers can be accommodated with modules available as either fully pregassed and tested systems or alternatively, mounting plates can be included for installation of the gas outlets by the medical gas contractor on site.

The units can be fed either via a bespoke surface mounted vertical riser section, through the back of the unit through the main body or through the side of the unit through the end cap.

Features

- These can be configured to incorporate any configuration of electrical sockets, medical gases, data and nurse call systems as well as any bespoke requirements.
- Modular design and assembly of the carcass of the unit that allows either the medical gas compartment or the data/power compartment to be increased in size. Where there is a space restriction or a large number of services, this can easily be accommodated.
- Fully configurable to suit every installation
- Dialysis system integration with numerous RO providers
- Bedhead units are supplied fully assembled, pre wired and gassed where required to ensure ease of installation
- Hygienic wipe down design thanks to clean, screw free fascia
- Trunking can be labelled and engraved to individual requirements
- Systems can be complemented with medical rail provide a complete solution



Healthcare Design.



Any type of nurse call can be easily integrated with provision for a bespoke nurse call plate which can also fit the required number of data outlets. The design of the units ensures correct separation is maintained between electrical, data and medical gases.



Professional Healthcare Design.

www.lipmedical.com

T: 01244 288 299
enquiries@lipmedical.com

Office

Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA

