

Mediwall, Vertical Bedhead Units



MANUFACTURED IN THE UK




LIP MEDICAL
HEALTHCARE DESIGN

Professional Healthcare Design.

www.lipmedical.com



T: 01244 288 299
enquiries@lipmedical.com

Our primary focus is the understanding of and commitment to enhance the patient environment.

With a broad range of knowledge and experience, we feel passionately about both design and functionality of the patient environment.

**We are
LiP Medical**

Healthcare Design.

Medilwall



Supply Units

Vertical Bedhead Supply Unit

Any nurse call supplier can be integrated within either horizontal or vertical units.

A full range of electrical sockets can be provided. Special locations such as A&E, Critical Care and Operating Departments can be provided with blue IPS sockets which are wired accordingly. Standard ward areas can be fitted with a range of white sockets with either white or red switches.

Medical gas terminal units from all manufacturers can either be provided as a fully pregassed and tested system or alternatively, mounting provision can be included for installation on site.


Both horizontal and vertical trunking can include horizontal medical rail. This can be supplied as part of a larger equipment rail system for areas such as intensive care, neonatal or A&E.

Anglepoise reading lights can be included on either horizontal or vertical units where required.

Vertical units are fed from the ceiling with cables and conduits routed down each side channel to ensure separation between mains and low voltage.

Accessories such as medical rail, infusion poles and examination/reading lighting brackets can be added as required and where a large number of services are required, a 'Wide' version is available to ensure that the units can meet all project requirements.





The Mediwall is a vertical bedhead unit which can be tailored to suit individual user requirements and offers a simple, stylish solution to deliver all the essential services to the patient. Units are supplied fully assembled and can either be pre wired for power or pre piped for medical gases off site to minimise installation time.

Healthcare Design.



Professional Healthcare Design.

www.lipmedical.com

T: 01244 288 299
enquiries@lipmedical.com

Office

Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA

