

**LP-MEDLN:**

# Mediline Horizontal Bedhead Unit

[www.lipmedical.com](http://www.lipmedical.com)



**MANUFACTURED IN THE UK**



  
**LIP MEDICAL**  
HEALTHCARE DESIGN

---

## Specifications:

# Mediline Horizontal Bedhead Unit

<b>Electrical Socket Type</b>	13A Switched/unswitched BS single or twin sockets - Quantity to client specification
<b>Equipotential Sockets</b>	POAG ID-6 as standard to client specification
<b>Data Outlets</b>	Single/Twin available to client specification
<b>Gas Outlets</b>	Quantities to suit project requirements - Pneumatech East/Zeus or Medaes gem 10 as standard
<b>Nurse Call Provision</b>	Compatible with all nurse call manufacturers
<b>Construction</b>	Powder coated aluminium extrusion with powder coated steel accessory plates to suit. Flexible to suit any number of beds with overall lengths designed to project requirements, supplied fully assembled with jointing kits where longer runs are required.
<b>Service Entry</b>	Back fed through carcass or side fed through end cap. A bespoke service riser is also available for surface service entry.
<b>Colour</b>	RAL9010 as standard. Optional colours available on request
<b>Assembly</b>	Supplied as fully assembled modules as standard with optional pre gassing and prewiring where required

---

T: 01244 288 299  
enquiries@lipmedical.com

### Office

Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA

