

LP-MEDLN:

Mediline Horizontal Bedhead Unit

www.lipmedical.com



MANUFACTURED IN THE UK




LIP MEDICAL
HEALTHCARE DESIGN

Specifications:

Mediline Horizontal Bedhead Unit

Electrical Socket Type	13A Switched/unswitched BS single or twin sockets - Quantity to client specification
Equipotential Sockets	POAG ID-6 as standard to client specification
Data Outlets	Single/Twin available to client specification
Gas Outlets	Quantities to suit project requirements - Pneumatech East/Zeus or Medaes gem 10 as standard
Nurse Call Provision	Compatible with all nurse call manufacturers
Construction	Powder coated aluminium extrusion with powder coated steel accessory plates to suit. Flexible to suit any number of beds with overall lengths designed to project requirements, supplied fully assembled with jointing kits where longer runs are required.
Service Entry	Back fed through carcass or side fed through end cap. A bespoke service riser is also available for surface service entry.
Colour	RAL9010 as standard. Optional colours available on request
Assembly	Supplied as fully assembled modules as standard with optional pre gassing and prewiring where required

T: 01244 288 299
enquiries@lipmedical.com

Office

Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA

