

**LP-MEDWLS:**

# Mediwall Vertical Bedhead Unit

www.lipmedical.com



MANUFACTURED IN THE UK



  
HEALTHCARE DESIGN

---

## Specifications:

# Mediwall Vertical Bedhead Unit

<b>Electrical Socket Type</b>	13A Switched/unswitched BS single or twin sockets - Quantity to client specification
<b>Equipotential Sockets</b>	POAG ID-6 as standard to client specification
<b>Data Outlets</b>	Single/Twin available to client specification
<b>Gas Outlets</b>	Quantities to suit project requirements - Penlon East/Zeus or Medaes gem 10 as standard
<b>Nurse Call Provision</b>	Compatible with all nurse call manufacturers
<b>Construction</b>	Powder coated aluminium extrusion with powder coated steel accessory plates to suit. Length to suit project requirements. Services entry from the top of the unit with gland points for electrical/data services and separate medical gas entry provided
<b>Colour</b>	RAL9010 as standard. Optional colours available on request
<b>Assembly</b>	Supplied as fully assembled modules as standard with optional pre gassing and prewiring if required

---

T: 01244 288 299  
enquiries@lipmedical.com

### Office

Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA

