
Case Study

Salford Royal Hospital



www.lipmedical.com

Professional Healthcare Design.

www.lipmedical.com



T: 0113 251 5119
F: 0113 251 5100
enquiries@lipmedical.com

Case Study

Salford Royal Hospital

Overview

Recognised as one of the best hospitals in the NHS, Salford Royal Hospital aims to deliver the safest and best patient care. Located in Manchester and with in excess of 7000 people working at the trust, Salford Royal is one of the larger hospitals within the Greater Manchester area.

As part of an ongoing refurbishment of their existing departments, LiP medical were invited to develop a solution for the bedhead supply units in the Operating Theatres recovery area during the latest phase of refurbishment.

Within the recovery area one of the most important factors to the clinical and design team was a space that remained functional but remained in keeping with the interior design to provide a less clinical feel.

The Solution

The Milan

The Milan unit from LiP Medical was chosen by both client team and design team as it offered the maximum flexibility to fit the design of the space whilst delivering the essential services to the patient. LiP Medical worked closely with the project team to design a solution that met and exceeded their requirements.

Healthcare Design.



The Milan



The Specifics

The Milan unit is part of the larger range of enhanced supply units that LiP Medical provide with the intention of improving the patient space.

The unit is designed around an aluminium framework that can be configured to suit any size and which supports the medical supply units.

The medical supply units are individually designed based on the project requirements and the unit is finished with HPL high pressure laminate panels which can be selected from a large range of finishes to compliment the interior design



Essential Services

Whilst aesthetics are a key factor in the design of the patient space, the fundamental reason for the bedhead supply unit cannot be forgotten especially in a recovery area where there is a much higher level of care required than in a normal ward. The units must always be designed to meet the safety standards whilst maintaining their functionality and overall ergonomics.

The Milan unit integrates two fully compliant medical supply units either side of the bed which are individually configured to house the power, data, nurse call and medical gases to the project teams requirements. The Milan unit at Salford Royal was configured specifically to include low level ventilation which is now a standard requirement in recovery areas.

The Accessories

A large number of accessories are available to compliment the Milan Panels.

As general lighting for the room and local illumination for the patient was required, the OKI-Fluo lamp was selected and fitted with matching laminate front strip to seamlessly integrate to the panel.

Bed buffers colours were selected to match the supply unit colour and fitted to the unit to prevent any risk of damage to the bedhead units.

Healthcare Design.





The Details

LiP Medicals broad portfolio of enhanced bedhead units allow us to design and offer a number of different solutions.

Working in conjunction with designers is key to our philosophy to achieve a warmer, less clinical environment for the patient.

Complimented by the large range of accessories that can be included with the bedhead units, this provides a comprehensive solution that can be adapted to suit each individual project.

Professional Healthcare Design.

www.lipmedical.com

T: 0113 251 5119

F: 0113 251 5100

enquiries@lipmedical.com

Offices

1200 Century Way, Thorpe Park Business Park, Colton, Leeds, LS15 8ZA

Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA