

LP-DIASYS:

'DIASYS' LED Minor Ops Lamp

www.lipmedical.com



The future looks good.


LIP MEDICAL
HEALTHCARE DESIGN

Specifications:

PR-DIASYS

Light Source: LED

Maximum Central Lighting*: 65,000 LUX

Adjustable Luminous Intensity*: From 30,000 LUX to 65,000 LUX

Fixed Field Diameter*: 14cm

Depth of Field L1+L2 (60% Ec) 2008*: 128cm

Colour Rendering Index Ra*: 96

Red Colour Rendering R9*: 95

Colour Temperature*: 4200 K

LED Life Time: 60,000 h

Max. Consumption: 7w

Dimension of the Light: 30 x 24.5cm



T: 01244 288 299
enquiries@lipmedical.com

Office

Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA

