

LP-AMLED:

'HUGO Ambience' LED Wall Mounted Over Bed Light



The future looks good.


LIP MEDICAL
HEALTHCARE DESIGN

Specifications:

LP-AMLED

Lamp Type	LED module, 12W
Light Intensity	50,000 Lux@0.5m
Light Field Diameter	Ø150mm@0.5m
Colour temperature	4,500 Kelvin
LED Life Expectancy	>45,000 hours
Control	Illuminated On/Off Switch mounted on the lamp head
Safety Standards	EN60601-1, EN60601-2-41
Arm Length	1540mm
Supply Voltage	230V, 50-60Hz
Total Weight	Approx. 7.5kg
Class of Protection	Class II
Floor to Ceiling Height	
Range:	2500mm - 2700mm 2700mm - 2900mm 2900mm - 3100mm

T: 01244 288 299
enquiries@lipmedical.com

Office

Unit 60 AML Zone 2, Deeside Industrial Estate, Third Avenue, Deeside, Flintshire, CH5 2LA

